

# mytardis

Connecting scientific instrument data to people

- Solving the problem of archival, archival and citation of raw data
- Easy sharing, access and publication of terabyte datasets
- More than 2 million files, 20 terabytes collected automatically for access - and growing
- 18 code contributors (python, [github.com/mytardis](https://github.com/mytardis))

MyTardis Home About Data Stats Register Log In


Experiment  
**X-ray crystal structure of the streptococcal specific phage lysin PlyC**  
Sheena McGowan , Ashley Buckle , James Whisstock

2 531 10.3 GB today Public

Hide Description

DOI for this data: [10.5072/03/50171CC7EEC1D](https://doi.org/10.5072/03/50171CC7EEC1D)

Sample diffraction image (PlyCB: plycb0001.img):



2 Datasets

Download Selected

PlyCB 120 3.8 GB

PlyC 411 6.5 GB

contact  
[steve.androulakis@gmail.com](mailto:steve.androulakis@gmail.com)

## Search

Experiment

## X-Ray Diffraction Data

Steve Androulakis

1

3

54.0 MB

28th September 2012

Public

## Describe

Toggle Full Description

A test of diffraction data such as that from the macromolecular beamlines at the Australian ...

Description

Metadata

Sharing

Transfer Datasets

## Share

Institution Monash University

Licensing This experiment data is licensed under  
Creative Commons Attribution 3.0  
Australia (CC BY 3.0).



Administrators Mister Mytardis

Download All

ZIP


TAR

## Publish

## 1 Dataset

Download Selected

Just start typing to filter datasets based on description



Diffraction Dataset

Download

3

54.0 MB

and Provide..

## Dataset

# Diffraction Dataset


From the experiment: [X-Ray Diffraction Data](#)

1

3

54.0 MB

## Preview Images



xtal2\_2\_4\_008.img

## 3 Files


[Download Selected Files](#)Select:  / 

<input type="checkbox"/>	<a href="#">xtal2_2_4_008.img</a> (18.0 MB)			
<input type="checkbox"/>	<a href="#">xtal2_2_4_006.img</a> (18.0 MB)			
<input type="checkbox"/>	<a href="#">xtal2_2_4_036.img</a> (18.0 MB)			



### TRD datafile

Preview Image



Detector Distance 299.947815 mm

Detector S/N 928

Direct Beam X-Pos 157.942017 mm

Direct Beam Y-Pos 160.654007 mm

Exposure Time 1.0 s

Image Size X 3072.0 px

Image Size Y 3072.0 px

Two Theta 0.0 °

## Public Access

### Step 1: Change Public Access:

Public access

Public

### Step 2: Select a license:

Use

**Creative Commons Attribution 3.0 Australia (CC BY 3.0)**

This licence lets others distribute, remix, tweak, and build upon your work, even commercially, as long as they credit you for the original creation. This is the most accommodating of licences offered under Creative Commons.



Selected

**Creative Commons Attribution-ShareAlike 3.0 Australia (CC BY-SA 3.0)**

This licence lets others remix, tweak, and build upon your work even for commercial purposes, as long as they credit you and licence their new creations under the identical terms.



## Sharing

### Users

Users who have a share in this experiment:

Username	Name	Permissions
mytardis	Steve Androulakis	<span>Read</span> <span>Edit</span> <span>Owner</span>
oded	Oded Kleifeld	<span>Read</span>
bosco	Bosco Ho	<span>Read</span>
synchrotron	Synchrotron Test	<span>Read</span> <span>Edit</span> <span>Owner</span>

[Change User Sharing](#)

### Groups

Groups who have a share in this experiment:

*There are currently no groups with access to this experiment.*

[Change Group Sharing](#)

## Links

This experiment is private. A temporary link can be created by its owner(s) and privately shared for direct access.

Temporary access links provide full access to recipients regardless of an experiment's public status.

Temporary Link	Expiry	Granted By
<a href="#">Right click.. copy link</a>	15th November 2012	mytardis
<a href="#">Right click.. copy link</a>	7th December 2012	mytardis

[Create New Temporary Link](#)

## User Sharing

<b>Oded Kleifeld [oded]</b> (oded.kleifeld@monash.edu)	<span>Read</span>	<span>Edit</span> <span>Delete</span>
<b>Bosco Ho [bosco]</b> (Bosco.Ho@monash.edu)	<span>Read</span>	<span>Edit</span> <span>Delete</span>
<b>Synchrotron Test [synchrotron]</b> (Tom.Caradoc-Davies@synchrotron.org.au)	<span>Read</span> <span>Edit</span> <span>Owner</span>	<span>Edit</span> <span>Delete</span>

### Add new user

Username:

Permissions:  
View Only

*Owners have the ability to change access controls and share experiments with others.*

[+ Add User](#)




# capturing data at its origin

Australian Synchrotron




MyTardis @ Monash University




Make Public


Cited In Papers



- Our scientists regularly generate large amounts of data at the Australian Synchrotron, so in collaboration, we connected their instruments to MyTardis for automatic capture of data.
- 4 nodes were established around Australia, servicing 8 universities. However, Monash University have recently made a deal with the Australian Synchrotron to store and mediate access to all users' data (starting with MX), regardless of their originating institution.

- 3rd party transfer between lab's network storage drive and MyTardis
- User begins transfer via MyTardis
- User receives email when transfer has completed and metadata extracted / processed





Users of the Hitachi H7500 TEM have their data automatically collected in the web by MyTardis.

MyTardis stores their data privately with sharing options for collaborators and eventual publication.

Data is available for download within the CVL desktop.



# Supporting Bioformats

<http://www.openmicroscopy.org/site/support/bio-formats/supported-formats.html>

<a href="#">Burleigh</a>	.img	■ ▼ ▼ ▼ ▼ ✕ ✕
<a href="#">Canon DNG</a>	.cr2, .crw	■ ■ ▼ ▼ ▼ ✕ ✕
<a href="#">Cellomics</a>	.c01	▲ ▼ ▼ ▼ ▼ ✕ ✕
<a href="#">cellSens VSI</a>	.vsi	▼ ■ ▼ ▼ ▼ ✕ ✕
<a href="#">DeltaVision</a>	.dv, .r3d	▲ ■ ■ ■ ■ ✕ ✕
<a href="#">DICOM</a>	.dcm, .dicom	▲ ▲ ▲ ■ ▼ ✕ ✓
<a href="#">ECAT7</a>	.v	■ ■ ▼ ▼ ▼ ✕ ✕
<a href="#">EPS (Encapsulated PostScript)</a>	.eps, .epsi, .ps	■ ■ ■ ▲ ▼ ✓ ✓
<a href="#">Evotec/PerkinElmer Opera Flex</a>	.flex, .mea, .res	▲ ▲ ▼ ▼ ▼ ✕ ✕
<a href="#">FEI</a>	.img	▼ ▼ ▼ ▼ ▼ ✕ ✕
<a href="#">FEI TIFF</a>	.tiff	▲ ■ ■ ▼ ▼ ✕ ✕
<a href="#">FITS (Flexible Image Transport System)</a>	.fits	▲ ▼ ▲ ■ ▼ ✕ ✓
<a href="#">Gatan Digital Micrograph</a>	.dm3	▲ ■ ▼ ▼ ▼ ✕ ✕
<a href="#">Gatan Digital Micrograph 2</a>	.dm2	■ ▼ ▼ ▼ ■ ✕ ✕
<a href="#">GIF (Graphics Interchange Format)</a>	.gif	▲ ▲ ▼ ▲ ▼ ✕ ✓
<a href="#">Hamamatsu Aquacosmos NAF</a>	.naf	■ ▼ ▼ ▼ ▼ ✕ ✕
<a href="#">Hamamatsu HIS</a>	.his	■ ▼ ▼ ▼ ▼ ✕ ✕
<a href="#">Hamamatsu ndpi</a>	.ndpi	▼ ■ ■ ▼ ▼ ✕ ✕
<a href="#">Hamamatsu VMS</a>	.vms	■ ■ ▼ ▼ ▼ ✕ ✕
<a href="#">Hitachi S-4800</a>	.txt, .tif, .bmp, .jpg	▲ ▲ ▲ ▼ ▼ ✕ ✕
<a href="#">ICS (Image Cytometry Standard)</a>	.ics, .ids	▲ ▲ ▲ ▲ ▲ ✓ ✓

## Datafile Metadata

For the file: test9.dm3 (2.2 MB)

### Bioformats Information

Preview Image






Image Attribute	Series #0 :
Image Attribute	Image count = 1
Image Attribute	RGB = false (1)
Image Attribute	Interleaved = false
Image Attribute	Indexed = false (true color)
Image Attribute	Width = 1024
Image Attribute	Height = 1024
Image Attribute	SizeZ = 1

# Growth




-  Synchrotron
-  Microscopy
-  Neutron Source
-  Medical Imaging






# Github: Contributions over time






# Github: Contributions over time



# Not so good

## no interactivity or social-ness



“Just let me get the files”



Steve  
Androulakis  
@bdp tardis

# Bioscience data platform

- A new, federally funded (NeCTAR) project involving MyTardis
  - Collaboration between Monash University and RMIT until 2014
  - NeCTAR is a federal government group establishing a Research Cloud
  - BDP Focuses on 3 areas:
  - Building a federated index of MyTardis data
  - Connecting data to supercomputing resources for processing
  - Providing interactive, rich publications online with the ability to manipulate data
- <https://www.nectar.org.au/bioscience-data-platform-tardis-cloud>



## Bioscience Data Platform: TARDIS in the Cloud

### Bioscience Data Platform: TARDIS in the Cloud

The Bioscience Data Platform (BDP) aims to bring existing computational systems together in a way that allows scientists to seamlessly work with data from capture through to publication. The platform will be backed by the Australian Characterisation Environment Virtual Lab to exploit the advantage of its uniquely powerful computational infrastructure. The platform will reach all areas of the structural biology workflow, from the inception of a research project through to scientific publication.

Typically, a scientist will collect data at an instrument facility (such as the Australian Synchrotron) to find their data catalogued and easily accessible on the storage cloud. The BDP will then provide a set of tools that link to high throughput computing infrastructure, potentially saving the scientist days of computational time and manual effort.

Rich visualisations of data stored in the cloud will be available, and peers will be able to annotate and draw attention to all aspects in a social, collaborative manner. These same tools will assist the peer review process by enabling researchers to share data with scientific journals and eventually the public in a secure and rich environment that tells the full story behind a discovery.

The Bioscience Data Platform will leverage the MyTardis data management system, hosted as a national service on cloud infrastructure, with linkages to major instrument facilities and high performance computing resources around Australia.

### Contact

Steve Androulakis <[steve.androulakis@monash.edu](mailto:steve.androulakis@monash.edu)>



# thanks

- Steve Androulakis
- [steve.androulakis@monash.edu](mailto:steve.androulakis@monash.edu)
- Resources
  - bioscience-data-platform.posterous.com
  - [github.com/bioscience-data-platform](https://github.com/bioscience-data-platform)
  - [github.com/mytardis/mytardis](https://github.com/mytardis/mytardis)
  - Twitter: @bdp\_tardis